

# Installation Data Sheet

## 128-Channel Motorized Commutator


### Features

- 128 signal channels, 16 reference channels, and +V, -V, and GND
- Maximum rotation speed is one revolution per second (RPS)
- Three speed ranges with variable speed and adjustable response
- Extremely-low actuation force: 1000  $\mu$ Nm
- Contactless actuation for motion sensor that is electrically and acoustically noiseless
- 120V or 220V operation


### Installation Schematic

Dimensions are in millimeters (mm)


**Note:** All dimensions are correct at the time of publication. Specifications are subject to change without notice. For planning pre-installation work, please contact Plexon to obtain the latest specifications.



**Electronic Drive Unit (EDU)**



**EDU Power Supply**











**Connector Pin Specifications for Harwin Connectors**

Male Connectors (Headstage Connector Housing – input):

GND = ground  
 NC = not connected  
 Ref = reference channel


Channel numbers shown are for example purposes only – color matching between top and bottom connectors is of primary importance

<p><b>Brown (1-15)</b></p> <p>A1 – Channel 1                  A3 – Channel 3                  A5 – Channel 5                  A7 – Channel 7                  A9 – (-V)                  A11 – GND                  A13 – NC                  A15 – Channel 9                  A17 – Channel 11                  A19 – Channel 13                  A21 – Channel 15                  A23 – (-V)                  A25 – GND</p>  <p>A2 – Channel 2                  A4 – Channel 4                  A6 – Channel 6                  A8 – Channel 8                  A10 – (+V)                  A12 – Ref A1                  A14 – NC                  A16 – Channel 10                  A18 – Channel 12                  A20 – Channel 14                  A22 – Channel 16                  A24 – (+V)                  A26 – Ref A2</p>	<p><b>Green (65-80)</b></p> <p>E1 – Channel 65                  E3 – Channel 67                  E5 – Channel 69                  E7 – Channel 71                  E9 – (-V)                  E11 – GND                  E13 – NC                  E15 – Channel 73                  E17 – Channel 75                  E19 – Channel 77                  E21 – Channel 79                  E23 – (-V)                  E25 – GND</p>  <p>E2 – Channel 66                  E4 – Channel 68                  E6 – Channel 70                  E8 – Channel 72                  E10 – (+V)                  E12 – Ref E1                  E14 – NC                  E16 – Channel 74                  E18 – Channel 76                  E20 – Channel 78                  E22 – Channel 80                  E24 – (+V)                  E26 – Ref E2</p>
<p><b>Red (16-32)</b></p> <p>B1 – Channel 17                  B3 – Channel 19                  B5 – Channel 21                  B7 – Channel 23                  B9 – (-V)                  B11 – GND                  B13 – NC                  B15 – Channel 25                  B17 – Channel 27                  B19 – Channel 29                  B21 – Channel 31                  B23 – (-V)                  B25 – GND</p>  <p>B2 – Channel 18                  B4 – Channel 20                  B6 – Channel 22                  B8 – Channel 24                  B10 – (+V)                  B12 – Ref B1                  B14 – NC                  B16 – Channel 26                  B18 – Channel 28                  B20 – Channel 30                  B22 – Channel 32                  B24 – (+V)                  B26 – Ref B2</p>	<p><b>Blue (81-96)</b></p> <p>F1 – Channel 81                  F3 – Channel 83                  F5 – Channel 85                  F7 – Channel 87                  F9 – (-V)                  F11 – GND                  F13 – NC                  F15 – Channel 89                  F17 – Channel 91                  F19 – Channel 93                  F21 – Channel 95                  F23 – (-V)                  F25 – GND</p>  <p>F2 – Channel 82                  F4 – Channel 84                  F6 – Channel 86                  F8 – Channel 88                  F10 – (+V)                  F12 – Ref F1                  F14 – NC                  F16 – Channel 90                  F18 – Channel 92                  F20 – Channel 94                  F22 – Channel 96                  F24 – (+V)                  F26 – Ref F2</p>
<p><b>Orange (33-48)</b></p> <p>C1 – Channel 33                  C3 – Channel 35                  C5 – Channel 37                  C7 – Channel 39                  C9 – (-V)                  C11 – GND                  C13 – NC                  C15 – Channel 41                  C17 – Channel 43                  C19 – Channel 45                  C21 – Channel 47                  C23 – (-V)                  C25 – GND</p>  <p>C2 – Channel 34                  C4 – Channel 36                  C6 – Channel 38                  C8 – Channel 40                  C10 – (+V)                  C12 – Ref C1                  C14 – NC                  C16 – Channel 42                  C18 – Channel 44                  C20 – Channel 46                  C22 – Channel 48                  C24 – (+V)                  C26 – Ref C2</p>	<p><b>Violet (97-112)</b></p> <p>G1 – Channel 97                  G3 – Channel 99                  G5 – Channel 101                  G7 – Channel 103                  G9 – (-V)                  G11 – GND                  G13 – NC                  G15 – Channel 105                  G17 – Channel 107                  G19 – Channel 109                  G21 – Channel 111                  G23 – (-V)                  G25 – GND</p>  <p>G2 – Channel 98                  G4 – Channel 100                  G6 – Channel 102                  G8 – Channel 104                  G10 – (+V)                  G12 – Ref G1                  G14 – NC                  G16 – Channel 106                  G18 – Channel 108                  G20 – Channel 110                  G22 – Channel 112                  G24 – (+V)                  G26 – Ref G2</p>
<p><b>Yellow (49-64)</b></p> <p>D1 – Channel 49                  D3 – Channel 51                  D5 – Channel 53                  D7 – Channel 55                  D9 – (+V)                  D11 – GND                  D13 – NC                  D15 – Channel 57                  D17 – Channel 59                  D19 – Channel 61                  D21 – Channel 63                  D23 – (-V)                  D25 – GND</p>  <p>D2 – Channel 50                  D4 – Channel 52                  D6 – Channel 54                  D8 – Channel 56                  D10 – (+V)                  D12 – Ref D1                  D14 – NC                  D16 – Channel 58                  D18 – Channel 60                  D20 – Channel 62                  D22 – Channel 64                  D24 – (+V)                  D26 – Ref D2</p>	<p><b>Gray (113-128)</b></p> <p>H1 – Channel 113                  H3 – Channel 115                  H5 – Channel 117                  H7 – Channel 119                  H9 – (-V)                  H11 – GND                  H13 – NC                  H15 – Channel 121                  H17 – Channel 123                  H19 – Channel 125                  H21 – Channel 127                  H23 – (-V)                  H25 – GND</p>  <p>H2 – Channel 114                  H4 – Channel 116                  H6 – Channel 118                  H8 – Channel 120                  H10 – (+V)                  H12 – Ref H1                  H14 – NC                  H16 – Channel 122                  H18 – Channel 124                  H20 – Channel 126                  H22 – Channel 128                  H24 – (+V)                  H26 – Ref H2</p>


**Female Connectors (Preamp Cables -- output):**

GND = ground  
 NC = not connected  
 Ref = reference channel


**Brown (1-16)**



**Green (65-80)**




**Red (16-32)**



**Blue (81-96)**



**Orange (33-48)**



**Violet (97-112)**



**Yellow (49-64)**



**Gray (113-128)**



## Installation

Installation of the Plexon 128-Channel Motorized Commutator consists of three major steps:

- 1 Mounting the commutator to the enclosure or experimental arena
- 2 Connecting the components
- 3 Adjusting the settings and testing the installation

### 1) Mounting the commutator plate

- 1 Carefully review the installation schematic on Page 1 to make sure there is adequate clearance for the commutator body, the wire bundles, and the torque-arm rotation. To minimize interference, ensure adequate clearance for the motor power cable and route it as far as practical from the signal cables.
- 2 In the experimental enclosure, provide an opening and mounting holes that match the dimensions in the installation schematic. Choose a mounting location that, when combined with the length of the headset cable, provides the animal with adequate range of movement for the requirements of the experiment.
- 3 Secure the commutator to the surface of the enclosure with suitable fasteners.
 

**Note:** The commutator uses a stepper motor to provide incremental movement. The vibration from the stepper motor can eventually loosen normal fasteners. To secure the commutator, consider the use of prevailing-torque fasteners (nylon-insert lock nuts, stover nuts, etc.) or a thread-locking compound such as Loctite®.
- 4 When mounting the commutator is completed, proceed with connecting the components.



### HINT


#### Avoiding electrical noise on signal channels

To avoid picking up unwanted ambient noise from the lab environment, the commutator must be properly grounded. Make sure to connect the green ground wire specified in [Step 1 on page 4](#) (in **2 Connecting the components** below) and use the terminating ground plugs for all unused channels as specified in [Step 6 on page 7](#).

### 2) Connecting the components

- 1 Attach the green ground wire between the ground screw on the commutator mounting base and the ground screw at the rear of the commutator control box. Make sure to use the

ground screw at the left rear corner of the box, which has the paint removed around it to provide a good ground.



- 2 Attach the AC cord to the power supply. At the **POWER** connector, make sure the power switch is **OFF**. Plug in the power supply. Connect the power supply output cable to the **POWER** connector on the rear panel of the EDU.






**CAUTION**


**Connector Lock Release**

The control box connectors for the SENSE IN and STEP MOTOR cables feature spring-loaded locks on each side. These locks secure the cable by grasping a raised tab on the cable connector. Before removing these cables from the control box, use a small blunt instrument such as a blade screwdriver to release the locks. Avoid damage to the connectors. Before removing the cables, always release the locking tabs.




**Connector lock release tab. Push out on each side to release cable.**

- 3 Ten external (outside) cables connect the commutator to the control box and the preamp: eight signal cables, a motor power cable, and a motion-sensor cable. At the back panel of the control box, make sure the beveled edges on the connector face up, then connect the motor power cable to the **STEP MOTOR P1** connector.



- 4 Make sure the beveled edges on the connector faces up and to the right, then connect the motion sensor cable to the **SENSE IN P2** connector.




- 5 Connect each color-coded data cable to its mating board in the preamp. The following table lists the cable colors for each channel range:

Color	Channel Range
Brown	1 to 16
Red	17 to 32
Orange	33 to 48
Yellow	49 to 64
Green	65 to 80
Blue	81 to 96
Violet	97 to 112
Gray	113 to 128

**Note:** When routing the cables, keep the data cables as far as possible from the motor power cable.

- 6 Eight internal (input) signal cables connect the headstages to the commutator. Inside the experimental enclosure, plug the headstage cables into the matching colored connectors on the commutator connector housing; see “Connector Pin Specifications” on page 2 The


following illustration shows the connector housing. (The connectors shown are female connectors, but the product “as shipped” has the connector gender changed to male.)



**NOTE:** The area within the red oval is not used in this commutator and is not described in this document.

The connectors are color-coded to match the previously described connectors/cables.

**Top Connectors**  
brown, orange, green (top to bottom)  
violet (right side)



**Bottom Connectors**  
blue, yellow, red (top to bottom)  
gray (left side)

**Area Within Green Oval**  
The connector is the stimulation pass-through connector  
The white identifiers at the left are:  
G = ground  
+ = +V  
- = -V

If all the cables are not connected to the connector housing, attach a terminating ground plug to each open connector. The terminating ground plug completes a path to ground for the “open” channels, which is essential to reduce noise on all signal channels.



- 7 Connect the other ends of the cables to the headstages.

Connection testing

- 8 Connect the headstages to Plexon headstage tester units. Connect the headstage tester units to the computer and provide a signal to the headstages.

**Note:** For information on the Plexon headstage tester unit, refer to the Plexon *RASPUTIN User's Guide*.

- 9 Turn the computer on and set up RASPUTIN Sort Client or Recorder to receive signals. On the MAP unit, turn on the preamp power switch. Verify that Sort Client or Recorder receives the test signals from the headstage tester unit. If Sort Client or Recorder do not show the signals, connect a headset cable directly to the preamp to verify that signals are being transmitted. Reconnect the headset cable to the commutator. If the signals still do not appear in Sort Client or Recorder, contact Plexon customer service for assistance.

- 10 Disconnect the headstage tester units.

End of connection testing

- 11 Turn the preamp power off.

3) Adjusting settings and testing the installation

- 1 On the EDU front panel, set the motor function switch to **DISABLE**.



- 2 At the **POWER** connector, turn the power **ON**.



- 3 On the front panel, set both the **MOTOR STEP SIZE** switches **MS1** and **MS2** to 1.

- 4 Set **MAX MOTOR SPEED** fully counter-clockwise to the slowest setting. Set **ADJ RESPONSE** fully counter-clockwise to the lowest setting.



### CAUTION

#### Possible Cable Damage

If the **SENSOR** cable is unplugged, the commutator turns independently and it can quickly twist the headset cables. Before setting the **MOTOR** switch to **ENABLE**, *always* make sure the headset cables are visible. If the commutator begins to twist the cables, immediately set the **MOTOR** switch to **DISABLE** and check the cable connections.

- 5 Set the **MOTOR** switch to **ENABLE**. Turn **ADJ RESPONSE** clockwise to increase the response at the connector housing. Gently rotate the housing back and forth to gauge the

response. Set the **ADJ RESPONSE** control to allow some connector housing movement in each direction without engaging the motor.

- 6 Turn the **MAX MOTOR SPEED** clockwise to increase the rotation speed. Gently rotate the connector housing toward its full travel in one direction to gauge the rotation speed. Set the **MAX MOTOR SPEED** to the desired setting for the experiment.
- 7 The rotation speed may also be changed by changing the step size. The various positions of the two **MOTOR STEP SIZE** switches provide three settings. Set both switches to **1** for the smallest step size, and hence the slowest and smoothest rotation. Set both switches to **0** for the largest step size, and hence the fastest rotation speed. Use **MAX MOTOR SPEED** to fine tune the speed in each range. See the following chart for the motor speed settings:

<b>Speed Range Switch Positions MS1 - MS2</b>	<b>Approximate Max Rotation Speed (rps)*</b>
1 - 1 (1/8 step)	1/4 rps
0 - 1 (1/4 step)	1/2 rps
1 - 0 (1/2 step)	1 rps
0 - 0 (1/2 step)	1 rps

**The optimum motor speed can:**

- Reduce acoustical noise produced by the motor shaking the experimental enclosure (i.e., metal cages)
- Reduce electrical noise produced by mechanical vibration of components
- More closely match animal behavior

\* varies somewhat depending on ADJ RESPONSE setting

- 8 To test the response and the rotation speed, attach the headstages to the experimental animal and observe the performance of the commutator. Adjust the response and rotation speed to provide the lowest level of interference for the animal.

**Note:** Slower speed ranges use smaller steps, which reduce electrical noise and mechanical vibration. In general, use the slowest speed range that is adequate for the experiment, then do a fine speed adjustment with the MAX MOTOR SPEED control.