

## **MEA Workstation Change Log (Last updated 2/12/07)**

### **Version 1.3.0**

#### **New Features:**

- Added support for Windows XP theme, minor display tweaks
- New Display Pause function to freeze the displays but still acquire data
- Added toolbar button for Assisted Sort view, added n, p, del keystroke handlers
- The licensing dialog makes it more obvious when multiple keys exist

#### **Improvements:**

- The Drop Count section of the status bar will flash red if the drop counts are non-zero
- Can now cancel out of an Add Unit in template sorting mode by hitting Esc
- If you add a new unit from the cluster displays by holding down Ctrl and left-clicking, the display will now flip into viewing the template set