

Alpha Omega Instruments Adaptors

Plexon-Alpha Omega Instruments Adaptors

Part: ADP/8o50-Alpha Omega-4-channel

Description: Adaptor to interface Plexon HST/8o50-GR headstage to 4 electrodes in an Alpha Omega drive using Multi-Contact LS1-B pins. 6 inch toothed clip on ground.



Provides an interface between: Plexon headstage

- HST/8o50-G1
- HST/8o50-G20

and Alpha Omega drive

- four electrodes in an Alpha Omega drive using Multi-Contact LS1-B pins. 6 inch toothed clip on ground

Not Pictured Part: ADP/8o50-AlphaOmega-8-channel

Description: Adaptor to interface Plexon HST/8o50-GR headstage to 8 electrodes in an Alpha Omega drive using Multi-Contact LS1-B pins. 6 inch toothed clip on ground.

Provides an interface between: Plexon headstages

- HST/8o50-G1
- HST/8o50-G20

and Alpha Omega drive

- eight electrodes in an Alpha Omega drive using Multi-Contact LS1-B pins. 6 inch toothed clip on ground