

Headstage and Gain	Channel #	Part #	Weight (grams)	Dimensions L x W x H (mm)	Picture and Technical Specifications
--------------------	-----------	--------	----------------	---------------------------	--------------------------------------

**Plexon  
4-  
channel  
EMG  
20x Gain**

4

HST/4EMG-  
8o50-G20

2.3

38.1 x 16 x  
3



**HST/4EMG-8o50-G20**


**Plexon  
12-  
channel  
EMG  
20x Gain**

12

HST/12EMG-  
26H-G20

5.9

30.5 x 29. x  
6



**HST/12EMG-26H-G20**