


Plexon Headstage Connector document



Headstage Type:
HST/8o50-G1-GR
HST/8o50-G1-TR
HST/8o50-G20-GR
HST/8o50-G20-TR

	Plexon Part Number	Omnetics Part Number
Mating Connector	CON/8o50m-10P	A11356-001
Input Connector	CON/8o50f-10P	A11862-001

Headstage Type:
HST/8o25-10P-GR
HST8o25-10P-TR

	Plexon Part Number	Omnetics Part Number
Mating Connector	CON/8o25m-10P	A8393-001
Input Connector	CON/8o25f-10P	A8376-001

Headstage Type:
HST/8o25-9P-GR
HST/8o25-9P-TR

	Plexon Part Number	Omnetics Part Number
Mating Connector	CON/8o25m-9P	A8177-001
Input Connector	CON/8o25f-9p	A7674-001

Headstage Type:
HST/16o25-18P-GR
HST/16o25-18P-TR

	Plexon Part Number	Omnetics Part Number
Mating Connector	CON/16o25m-18P	A8141-001
Input Connector	CON/16o25f-18P	A8104-001

Headstage Type:
HST/32o25-36P-GR
HST/32o25-36P-TR

	Plexon Part Number	Omnetics Part Number
Mating Connector	CON/32o25m-36P	A8648-001
Input Connector	CON/32o25f-18P	A8649-001

Headstage Type:
HST/16V-G20

	Plexon Part Number	Omnetics Part Number
Mating Connector	CON/16m-V	A8784-001
Input Connector	CON/16f-V	A8858-001

Headstage Type:
HST/32V-G20

	Plexon Part Number	Omnetics Part Number
Mating Connector	CON/32m-V	A8828-001
Input Connector	CON/32f-V	A8857-001