


## Adaptors

### Plexon-NAN Instruments Adaptors

#### Part: ADP/8o50-NAN-4

**Description:**

“Adaptor to interface 8-channel (10-pin) Omnetics 0.050 headstage to 1-4 NAN electrode drive channels; NAN tower mating pins on all channels, clip on ground, 4" channel and ground leads”



**Provides an interface between:** Plexon headstages

- HST/8o50-G1
- HST/8o50-G20


and NAN drive

- one to four electrodes

#### Part: ADP/8o50-NAN-8

**Description:**

Adaptor to interface 8-channel (10-pin) Omnetics 0.050 headstage to 5-8 NAN electrode drive channels.



**Provides an interface between:** Plexon headstage

- HST/8o50-G1
- HST/8o50-G20

and NAN drive

- five to eight electrodes