

Adaptors

Plexon-Tetrode Adaptors

Part:ADP/24TET-32V


Description: Chronic electrode adaptor for interfacing three 32-channel VLSI headstages to 24 tetrode wires; 6 pass throughlines for stimulation.

Provides an interface for: three Plexon headstages

- HST/32V-G20

to

- 24 tetrode wires
- 96 individual wires
- Neuralynx 12-Drive-H (Twelve Tetrode Microdrive)
- Kubie Drive



Not Pictured Part: ADP/32Vx4-NLynx12DriveH

Description: Adaptor to interface up to four 32-channel VLSI headstages to a Neuralynx 12-Drive-H.

Provides an interface for: four Plexon headstages

- HST/32V-G20

to

- 32 tetrode wires or 128 individual wires
- Neuralynx 12-Drive-H (Twelve Tetrode Microdrive)