

Adaptors

Plexon-Thomas Recording Adaptors

Part: ADP/TREC-5-PBX16-DB9f

Description: Adaptor cable to connect 5-channel Thomas Recording drive (with female DB9) to PBX/16 (male Harwin), 36" length

To Plexon Preamp -|- To 5-channel Thomas drive



Not Pictured Parts (cables)

- **ADP/TREC3-PBX16** - Adaptor cable to connect 3-channel Thomas Recording drive (with female DB25) to PBX/16 (male Harwin), 36" length
- **ADP/TREC5-PBX16-DB25f** - Adaptor cable to connect 5-channel Thomas Recording drive (with female DB25) to PBX/16 (male Harwin), 36" length
- **ADP/TREC7-PBX16** - Adaptor cable to connect 7-channel Thomas Recording drive (with female DB25) to PBX/16 (male Harwin), 36" length
- **ADP/TREC12-PBX16** - Adaptor cable to connect 12-channel Thomas Recording drive (with female DB25) to PBX/16 (male Harwin), 36" length
- **ADP/TREC16-PBX16** - Adaptor cable to connect 16-channel Thomas Recording drive (with female DB25) to PBX/16 (male Harwin), 36" length